# Earth field NMR

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Inductance  $\approx 73$  mH; d.c. Resistance  $\approx 10.5 \Omega$

Coil  $Q$ -factor (at 2.0 kHz):  $\approx 70$

Current Limit: 0.5 to 3.0 Amperes; 0.1 to 99.9 s duration

Notice that an *external* current-regulated power supply is needed for the sample coil used in its polarizing-coil mode.

##### **Sample holders:**

Both fit into the sample coil's tunnel of inside diameter 2.00"

Standard holder: 125 ml Nalgene bottle

Segmented holder: 7 compartments, 11 ml each

##### **Accessible frequencies:**

1.6 to 2.6 kHz

##### **Controller:**

Pre-amplifier gain:  $\approx 1000$

Band-pass amplifier gain:  $\approx 15$

### **Helmholtz coils:**

- 2 Coils, each of 30 turns of #20 AWG copper wire, connected in series
  - Nominal Diameter of coils: 604 mm
  - Nominal Separation of coils: 303 mm
  - Coil Constant: near 89 ( $\pm 1\%$ )  $\mu\text{T/A}$
  - Current Limit: 3.0 Amperes continuous current
  - Field Uniformity: 0.01 %, for volume within 3 cm of center of coils
  - Frame, permitting altitude/azimuth adjustment of coils' z-axis

Notice that an *external* current-regulated power supply is needed for the Helmholtz coils.

### **Gradient coils:**

- 3 Coil Systems, for creating gradients  $\partial B_z/\partial x$ ,  $\partial B_z/\partial y$ ,  $\partial B_z/\partial z$ 
  - Gradient coil systems are embedded in the planes of the Helmholtz-coil forms.
  - Coil Constants: near 250  $\mu\text{T/m}$  per Ampere of coil current
  - Current Supply: 3 supplies, each  $\pm 20$  mA, with 10-turn manual controls
  - Gradient Capability:  $\pm 5$   $\mu\text{T/m}$ , for  $\pm 5$  turns from center position

Notice that the EF-NMR G/FC Controller *contains* the current-regulated power supplies for these coils.

### **Current-Monitor points:**

- Monitor Resistor in series with Helmholtz coil: 0.1  $\Omega$ , 1% tolerance, 25 W
- Monitor Resistors in series with each gradient coil: 100  $\Omega$ , 1% tolerance,  $\frac{1}{4}$  W

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
Figure 1-13 User Interface

### 4. Sample Rate ; Memory Depth

Displays the current sample rate and memory depth of the analog channel.

The sample rate and the memory depth will change along with the horizontal time base.

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What happens when a bottle of water is placed in a magnetic field  $B$ ? The answer differs in the short term and the long term. The 'long term' here means after several seconds. After enough time, the protons do tend to align with the field. But the degree of alignment is quite small, due to the weakness of the magnetic interaction. Given a proton gyromagnetic ratio of order  $100 \times 10^6 \text{ C/kg}$ , and a proton intrinsic (or spin) angular momentum of order  $\hbar \approx 10^{-34} \text{ J s}$ , we expect a proton magnetic moment of

$$\mu = \gamma L \approx (10^8 \text{ C/kg})(10^{-34} \text{ J s}) = 10^{-26} \text{ J/T} = 10^{-26} \text{ A} \cdot \text{m}^2.$$

14. Confirm that the units are correct in this equation.

Then even in a field of  $B = 1 \text{ T}$  (the field of a strong permanent magnet), the magnetic energy of interaction of a proton is limited to  $|U_{\text{mag}}| = \mu B \approx 10^{-26} \text{ J}$ .

This is *much* smaller than the room-temperature thermal 'energy of dis-orientation', which is of order  $k_B T$ , where  $k_B$  is Boltzmann's constant.

15. How big is  $k_B T$ , in Joules, at room temperature?

Because the magnetic-field interaction is so very weak, the protons' magnetic moments line up with the field only to a very small degree. To use quantum-mechanical terms, instead of a 50:50 mix of 'spin-up' and 'spin-down' protons, we expect an eventual equilibrium of  $50 + \epsilon$  :  $50 - \epsilon$  for spin-up vs. spin-down protons. But the difference in number of spin-up and spin-down protons leaves a (small) un-cancelled, or net, magnetic moment of the sample.

The way to quantify this is to compute the magnetization  $M$  of the sample, the net magnetic moment (in  $\text{A} \cdot \text{m}^2$ ) per unit volume (in  $\text{m}^3$ ). Thus  $M$  has the units of  $\text{A/m}$ . The value expected in our application is derived in Appendix A1, and gives

$$M = \frac{N}{V} \left( \frac{\gamma \hbar}{2} \right)^2 \frac{B}{k_B T}$$

Here  $N$  is the number of polarizable protons, and  $V$  is the volume they're contained in.

16. Compute the value of  $M$  that would result for water at room temperature, in a field of 50 mT.
17. Compute, for comparison, the value of  $M$  that would result if you could *fully align* the magnetic moments of all the protons in a water sample. Compare both these values to the number  $M \approx 1. \times 10^6$  A/m that applies to NdFeB rare-earth permanent-magnet material.